



REF CA-09803 8.5^{Fr.}



PSI for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue
- Dilator Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow 1: Advancer
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle 1:
- Pressure Transduction Probe 1:
- Injection Needle: 25 Ga. x 1st (2.54 cm) and 3 mL Luer-Lock Syringe Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® 1:
- Raulerson Spring-Wire Introduction Syringe Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip
- 1: Syringe Obturator: 8 Fr. 1:
- Cath-Gard[®] Catheter Contamination Shield: 80 cm with TwistLock™ Adapter 1:

- Adapter
 SharpsAway® Disposal Cup
 Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
 Drape: 24" x 36" (61 cm x 91 cm)

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. Caution: This Product Contains Dry Natural Rubber.

Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic.

- Scalpel: #11 1:
- 1: 2:
- Stopcock: 4-way High-Flow Gauze Pad: 2" x 2" (5 cm x 5 cm) Gauze Pad: 4" x 4" (10 cm x 10 cm) 5:
- 1: Suture: 3-0 Braided Silk with Straight Needle

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL058306 R00 (2021-11)



LOT

ςĮ

PSI for use with 7 - 7.5 Fr. Catheters

8.5^{Fr.}

