



REF AK-09903-A

9^{Fr.} 10^{cm sheath length} .035^{inch dia. spring-wire guide}

PSI for use with 7.5 - 8 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®¹ One-Step Solution with Hi-Lite Orange™ Tint
- 1: Obturator: 8 Fr.
- 1: Suture: 3-0 Braided Silk with Straight Needle
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 1: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Scalpel: #11
- 1: Stopcock: 4-way High-Flow
- 7: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Tubing±: 20" (50 cm) Extension

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov

±Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F). Avoid freezing and excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL058309 R00 (2022-03)



REF AK-09903-A

9^{Fr.} 10^{cm sheath length} .035^{inch dia. spring-wire guide}

PSI for use with 7.5 - 8 Fr. Catheters