



10 cm sheath length .035 inch dia. spring-wire guide

## PSI Kit for use with 7.5 Fr. Catheters

## Contents:

- Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer

- Other) with Arrow Advancer
  Catheter: 18 Ga. x 1-3/4" (4.45 cm) Radiopaque over 20 Ga. TW Introducer Needle
  Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 3 mL Luer-Lock Syringe
  Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
  Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Luer-Slip Syringe
  Obtractor 8 Er.
- Obturator: 8 Fr. Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter

- Cath-Garde Catheter Contamination Shield: 80 cm with 1 wistLock." Adapter Dust Cap; Non-Vented SharpsAway® Disposal Cup Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive Scalpel: #11 Stopcock: 4-way High-Flow Gauze Pad; 2" x 2" (5 cm x 5 cm)

- Gauze Pad: 2" x 2" (5 cm x 5 cm)
  Gauze Pad: 4" x 4" (10 cm x 10 cm)
  Suture: 3-0 Braided Silk with Straight Needle

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. Caution: This Product Contains Dry Natural Rubber.
Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.











Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL058314 R01 (2022-02)





