



8.5^{Fr.}

10 cm
sheath
length

.035 inch dia.
spring-wire
guide

Contents:

- | | |
|--------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------|
| 1: Arrow® and Blue® Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator | 1: SharpsAway® Disposal Cup |
| 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: Drape: 24" x 36" (61 cm x 91 cm) |
| 1: Injection Needle: 25 Ga. x 1" (2.54 cm) | 1: Drape: 41" x 55" with 4" fenestration |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: Needle Holder |
| 1: Pressure Transduction Probe | 1: Safety Scalpel: #11 |
| 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Stopcock: 4-way High-Flow |
| 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Syringe: 3 mL Luer-Lock | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | 1: Suture: 3-0 Silk with Curved Needle |
| 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep® ¹ One-Step Solution with Hi-Lite Orange™ Tint | 1: Suture: 3-0 Braided Silk with Straight Needle |
| 1: Obturator: 8 Fr. | 1: Tubing: 20" (50 cm) Extension |
| 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter | |
| 1: Dust Cap: Non-Vented | |
| 1: SharpsAway® II Locking Disposal Cup | |
- ¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

California Prop. 65
⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

±Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F)

Sterilized using ethylene oxide unless otherwise indicated in the contents list

Fluid path components are non-pyrogenic.

Contraindications: The Arrow®ard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).



LOT



 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LB 058321 R01 (2022-02)



8.5^{Fr.}

10 cm sheath length

.035 inch dia.
spring-wire
guide

