



REF AK-09903-CDC

9 Fr.

10 cm  
sheath  
length

.035 inch dia.  
spring-wire  
guide

# PSI for use with 7.5 - 8 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator

1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle

1: Injection Needle: 25 Ga. x 1" (2.54 cm)

1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

1: Pressure Transduction Probe

1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe

1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe

1: Syringe: 3 mL Luer-Lock

1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep

1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>1</sup> One-Step Solution with Hi-Lite Orange™ Tint

1: Obturator: 8 Fr.

1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter

1: Dust Cap: Non-Vented

1: SharpsAway® II Locking Disposal Cup

1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration

1: Drape: 24" x 36" (61 cm x 91 cm)

1: Drape: 41" x 55" with 4" fenestration

1: Needle Holder

1: Safety Scalpel: #11

1: Stopcock: 4-way High-Flow

2: Gauze Pad: 2" x 2" (5 cm x 5 cm)

5: Gauze Pad: 4" x 4" (10 cm x 10 cm)

1: Suture: 3-0 Silk with Curved Needle

1: Suture: 3-0 Braided Silk with Straight Needle

1: Tubing±: 20" (50 cm) Extension
- <sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only**  
**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.  
**California Prop. 65**  
⚠ **WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)  
**±Precaution:** Contains Phthalate: DEHP.  
**Warning:** The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.  
**Not made with natural rubber latex.**  
Store between 20 – 25°C (68 – 77°F).  
Fluid path components are non-pyrogenic.



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA  
**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)  
LBL058331 R01 (2022-02)



REF AK-09903-CDC

9 Fr.

10 cm  
sheath  
length

.035 inch dia.  
spring-wire  
guide

# PSI for use with 7.5 - 8 Fr. Catheters