



## **REF AK-09803-CDC**

8.5 Fr.

10 cm sheath length .035 inch dia. spring-wire guide

1: Needle Holder
1: Safety Scalpel: #11
1: Stopcock: 4-way High-Flow
2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
1: Suture: 3-0 Silk with Curved Needle
1: Suture: 3-0 Braided Silk with Straight Needle
1: Tubing±: 20" (50 cm) Extension
1: Tubing±: 20" (50 cm) Extension
1: A registered trademark of Care Fusion or one of it

Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration Drape: 24" x 36" (61 cm x 91 cm) Drape: 41" x 55" with 4" fenestration Needle Holder

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

## PSI for use with 7.5 Fr. Catheters

- Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle

- Needle

  Injection Needle: 25 Ga. x 1" (2.54 cm)

  Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

  Pressure Transduction Probe

  Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow®
  Raulerson Spring-Wire Introduction Syringe

  Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe

  Syringe: 3 mL Luer-Lock

  5 mL 1% Lidocaine HCI Solution and Alcohol Prep

  3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint

  Obturator: 8 Fr.

  Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter

  Dust Cap: Non-Vented

  SharpsAway® II Locking Disposal Cup

  SharpsAway® Disposal Cup

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F). Fluid path components are non-pyroge









Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL058355 R01 (2022-02)





10 cm sheath length .035 inch dia. spring-wire guide

