



REF ASK-09903-SM



9^{Fr.}

10^{cm}
sheath
length

.035^{inch dia.}
spring-wire
guide

PSI Kit for use with 7.5 - 8 Fr. Catheters

Contents:

- | | |
|--|--|
| 1: Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator | 1: Needle Holder |
| 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: Forceps |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: Safety Scalpel: #11 |
| 1: Pressure Transduction Probe | 1: Stopcock: 4-way High-Flow |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Checklist/CLIP Sheet |
| 1: Injection Needle: SafetyGlide™ ¹ 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe | 1: Sterile Procedure Sign |
| 1: Injection Needle: SafetyGlide™ ¹ 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Syringe: 10 mL Luer-Lock | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Obturator: 8 Fr. | 1: Dressing: Tegaderm® ² 10 cm x 12 cm |
| 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter | 1: Suture: 3-0 Braided Silk with Straight Needle |
| 1: Dust Cap: Non-Vented | 1: Suture: 3-0 Silk with Curved Needle |
| 1: SharpsAway® II Locking Disposal Cup | 1: Tubing±: 20" (50 cm) Extension |
| 1: SharpsAway® Disposal Cup | 1: HemoHopper® Fluid Receptacle |
| 1: Maximal Barrier Drape™ with 4" fenestration | |
| 1: Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive | ¹ A trademark of Becton, Dickinson and Company. |
| 1: Towel | ² A registered trademark of 3M Company. |

±Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.
Not made with natural rubber latex.

