



REF ASK-29903-UHN

(IPNXXXXXX)

ASK ARROW®

9 Fr.

10 cm sheath length

.035 inch dia. spring-wire guide

Arrowgard Blue® PSI Kit for use with 7.5 - 8 Fr. Catheters

**Contents:**

- | | |
|--|--|
| 1: Arrowgard Blue® Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator | 1: Towel |
| 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: Needle Holder |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: Safety Scalpel: #11 |
| 1: Pressure Transduction Probe | 1: Stopcock: 4-way High-Flow |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Checklist/CLIP Sheet |
| 1: Injection Needle: SafetyGlide™ ¹ 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe | 1: Sterile Procedure Sign |
| 1: Injection Needle: SafetyGlide™ ¹ 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Syringe: 10 mL Luer-Lock | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Obturator: 8 Fr. | 1: Dressing: Tegaderm® ² 10 cm x 12 cm |
| 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter | 1: Suture: 3-0 Silk with Curved Needle |
| 1: Dust Cap: Non-Vented | 1: Tubing±: 20" (50 cm) Extension |
| 1: SharpsAway® II Locking Disposal Cup | 1: HemoHopper® Fluid Receptacle |
| 1: SharpsAway® Disposal Cup | ¹ A trademark of Becton, Dickinson and Company. |
| 1: Maximal Barrier Drape™ with 4" fenestration | ² A registered trademark of 3M Company. |
| 1: Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive | |

±Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Contraindications: The Arrowgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

