



REF CDC-29803-XCN1A

8.5 Fr.

10 cm sheath length

.035 inch dia. spring-wire guide

ARROW®

ErgoPack<sup>®</sup>  
Complete

# Arrowg+ard Blue<sup>®</sup> PSI for use with 7.5 - 8 Fr. Catheters

## Contents:

- |   |   |
|---|---|
| 1: Arrowg+ard Blue <sup>®</sup> Arrow-Flex <sup>®</sup> Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator  | 1: Dust Cap: Non-Vented   |
| 1: Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow <sup>®</sup> GlideWheel <sup>™</sup> Wire Advancer | 1: SharpsAway <sup>®</sup> II Locking Disposal Cup                        |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle  | 1: SharpsAway <sup>®</sup> Disposal Cup                                   |
| 1: Injection Needle: SafetyGlide <sup>™</sup> 1 25 Ga. x 1" (2.54 cm)   | 1: Maximal Barrier Drape <sup>™</sup> with 4" fenestration                |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)   | 1: Towel  |
| 1: Pressure Transduction Probe  | 1: Needle Holder  |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow <sup>®</sup> Raulerson Spring-Wire Introduction Syringe  | 1: Safety Scalpel: #11  |
| 1: Injection Needle: SafetyGlide <sup>™</sup> 1 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe  | 1: Stopcock: 4-way High-Flow  |
| 1: Syringe: 3 mL Luer-Lock  | 1: Checklist/CLIP Sheet   |
| 2: 10 mL Luer-Lock Pre-Filled Saline Syringe  | 1: Sterile Procedure Sign   |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep  | 1: Medication Label: 1% Lidocaine   |
| 2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep <sup>®</sup> One-Step Solution with Hi-Lite Orange <sup>™</sup> Tint   | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)                                       |
| 1: Obturator: 8 Fr.   | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)                                     |
| 1: Cath-Gard <sup>®</sup> Catheter Contamination Shield: 80 cm with TwistLock <sup>™</sup> Adapter  | 1: Surgical Apparel: Impervious Gown                                      |
|   | 1: Dressing: Tegaderm <sup>®</sup> CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm) |
|   | 1: Transducer Cover: 14 cm x 147 cm                                       |
|   | 1: Surgical Apparel: Mask with Eye Shield                                 |
|   | 1: Surgical Apparel: Bouffant Cap   |
|   | 1: Suture: 3-0 Silk with Curved Needle                                    |
|   | 1: Tubing <sup>±</sup> : 20" (50 cm) Extension                            |
|   | 1: HemoHopper <sup>®</sup> Fluid Receptacle                               |

<sup>1</sup>A trademark of Becton, Dickinson and Company.

<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

<sup>3</sup>A registered trademark of 3M Company.

**±Precaution: Contains Phthalate: DEHP.**

**Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.**

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

**Contraindications: The Arrowg+ard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).**



LBL058481 R01 (2023-12)