



REF **CU-22122-F** (IPNXXXXXX)

12 Fr. **2** Lumen **16** cm catheter length **.035** inch dia. spring-wire guide

You-Bend™ Arrowgard Blue® Two-Lumen Hemodialysis Catheter for High Volume Infusions

Contents:

- 1: Two-Lumen Catheter: 12 Fr. (4.06 mm OD) x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowgard® Antimicrobial Surface Treatment, Extension Line Clamps origin:xxxxxxxxxxxxxxxxxxxx
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer origin:xxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle origin:xxxxxxxxxxxxxxxxxxxx
- 1: Pressure Transduction Probe origin:xxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe origin:xxxxxxxxxxxxxxxxxxxx
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm origin:xxxxxxxxxxxxxxxxxxxx
- 2: Dust Cap: Non-Vented origin:xxxxxxxxxxxxxxxxxxxx
- 1: Drape: 24" x 36" (61 cm x 91 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxx

All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.
 www.P65Warnings.ca.gov

‡MR Conditional

Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels. The Arrowgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs.



Lumen	Priming Volume* (mL)	Flow Rate (mL/min)	Venous Pressure Setting†† (mmHg)
Distal (12 Ga.)	1.6	220 370	100 200
Proximal (12 Ga.)	1.5	220 370	100 200

* Priming volumes are approximate and are done without accessories.
 †† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.



EC REP Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
CE 2797

All components are CE 2797 unless otherwise noted.

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

Set Product of xxx
 Packaged in xxx

Teleflex®
 www.teleflex.com

LBL058901 R02 (2024-05)



You-Bend™ Arrowgard Blue® Two-Lumen Hemodialysis Catheter for High Volume Infusions



REF **CU-22122-F**

12 Fr. **2** Lumen **16** cm catheter length **.035** inch dia. spring-wire guide