

ARROW



REF PS-01672

Two-Lumen PICC

Contents:

- 1: Two Lumen Catheter: 16 Ga. (1.76 mm OD) x 70 cm, Blue FlexTip®, Contamination Guard
- 1: Peel-Away Catheter: 14 Ga. x 2-1/2" (6.35 cm) Radiopaque over 15 Ga. TW Introducer Needle and 3 mL Luer-Slip Syringe
- 1: Skin Protectant Prep
- 2: Dust Cap: Non-Vented
- 2: Clamp: Extension Line
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Patient Information Booklet
- 1: Patient ID Card
- 1: Dressing: STATLOCK®¹ Catheter Stabilization Device

¹A registered trademark of C. R. Bard, Inc.

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.

Contraindications: The PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.

ARROW



REF PS-01672

Two-Lumen PICC

Contents:

- 1: Two Lumen Catheter: 16 Ga. (1.76 mm OD) x 70 cm, Blue FlexTip®, Contamination Guard
- 1: Peel-Away Catheter: 14 Ga. x 2-1/2" (6.35 cm) Radiopaque over 15 Ga. TW Introducer Needle and 3 mL Luer-Slip Syringe
- 1: Skin Protectant Prep
- 2: Dust Cap: Non-Vented
- 2: Clamp: Extension Line
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Patient Information Booklet
- 1: Patient ID Card
- 1: Dressing: STATLOCK®¹ Catheter Stabilization Device

¹A registered trademark of C. R. Bard, Inc.

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.

Contraindications: The PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.

| | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) | Pump Flow Rate‡ (mL/hr) |
|-------------------|----------------------|----------------------------|-------------------------|
| Distal (18 Ga.) | 0.51 | 323 | 1410 |
| Proximal (20 Ga.) | 0.48 | 125 | 2657 |

* Priming volumes are approximate and are done without accessories.
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.
 ‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

en sl es



EU Authorized Representative and Importer:

EC REP Teleflex Medical IDA Business and Technology Park Dublin Road, Athlone 2797 Co. Westmeath, Ireland

LOT



Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

Teleflex® www.teleflex.com

LBL059139 R01 (2021-08)

| | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) | Pump Flow Rate‡ (mL/hr) |
|-------------------|----------------------|----------------------------|-------------------------|
| Distal (18 Ga.) | 0.51 | 323 | 1410 |
| Proximal (20 Ga.) | 0.48 | 125 | 2657 |

* Priming volumes are approximate and are done without accessories.
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.
 ‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

en sl es



EU Authorized Representative and Importer:

EC REP Teleflex Medical IDA Business and Technology Park Dublin Road, Athlone 2797 Co. Westmeath, Ireland

LOT



Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

Teleflex® www.teleflex.com

LBL059139 R01 (2021-08)



16 Ga. 2 Lumen

70 cm catheter length

REF PS-01672 Two-Lumen PICC



16 Ga. 2 Lumen

70 cm catheter length

REF PS-01672 Two-Lumen PICC