<sup>4</sup>A registered trademark of Johnson & Johnson Corporation.





REF NA-04020-S1ASP (IPNXXXXXX)

20 Ga. 4.45 cm catheter length of the spring-wire guide

## Radial Artery Catheterization with Integral Needle Protection

- <sup>1</sup>A trademark of Becton, Dickinson and Company. <sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries. <sup>3</sup>A registered trademark of TIDI.

Rx only
Warning: Read all package insert warnings, precautions, and
instructions prior to use. Failure to do so may result in severe patient
injury or death. www.teleflex.com/IFU
California Prop. 65

\( \triangle \) WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

\( \triangle \) Precaution: Contains Phthalate: DEHP.
Warning: The results of certain animal experiments have shown
phthalates to be potentially toxic to reproduction. Proceeding from
the present state of scientific knowledge, risks for male premature
infants cannot be excluded in the case of long-term exposure or
application. Medical products containing phthalates should be used
only temporarily with pregnant women, nursing mothers, babies and
infants.
Not made with natural rubber latex.

Store between 20 - 25°C (68 - 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.
































**LOT** LotNumber







(01)XXXXXXXXXXXXXXX (17)YYMMDD (11)YYMMDD

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

**Teleflex**®

LBL059315 R02 (2024-05)



**REF NA-04020-S1ASP** 

20 Ga. 4.45 cm catheter length of the control of th



Radial Artery Catheterization with Integral **Needle Protection** 

