



REF NA-04020-X1A

(IPNXXXXXX)

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guide

# Arterial Catheterization

### Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide origin:xxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: SafetyGlide™ 1 25 Ga. x 1" (2.54 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Syringe: 3 mL Luer-Lock origin:xxxxxxxxxxxxxxxxxxxx
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep origin:xxxxxxxxxxxxxxxxxxxx
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint origin:xxxxxxxxxxxxxxxxxxxx
- 1: Wing Clip: Suture origin:xxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® Disposal Cup origin:xxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® II Locking Disposal Cup origin:xxxxxxxxxxxxxxxxxxxx
- 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration origin:xxxxxxxxxxxxxxxxxxxx
- 1: Towel: 13" x 18" (33 cm x 46 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Needle Holder origin:xxxxxxxxxxxxxxxxxxxx
- 1: Safety Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxx
- 1: Checklist Sheet origin:xxxxxxxxxxxxxxxxxxxx
- 1: Medication Label: 1% Lidocaine origin:xxxxxxxxxxxxxxxxxxxx
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Dressing: Tegaderm®<sup>3</sup> 10 cm x 12 cm origin:xxxxxxxxxxxxxxxxxxxx
- 1: Dressing: BIOPATCH®<sup>4</sup> origin:xxxxxxxxxxxxxxxxxxxx
- 1: Surgical Apparel: Bouffant Cap origin:xxxxxxxxxxxxxxxxxxxx
- 1: Surgical Apparel: Mask with Eye Shield origin:xxxxxxxxxxxxxxxxxxxx
- 1: Suture: 3-0 Silk with Curved Needle origin:xxxxxxxxxxxxxxxxxxxx
- 1: Tubing±: 6" Extension with One-Way Stopcock origin:xxxxxxxxxxxxxxxxxxxx

<sup>1</sup>A trademark of Becton, Dickinson and Company.  
<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.  
<sup>3</sup>A registered trademark of 3M Company.  
<sup>4</sup>A registered trademark of Johnson & Johnson Corporation.

### Rx only

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**

### California Prop. 65

**⚠ WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

### ±Precaution: Contains Phthalate: DEHP.

**Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.**

### Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



**LOT** Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

Set Product of xxx  
Packaged in xxx

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA



www.teleflex.com

LBL059316 R02 (2024-05)

REF NA-04020-X1A

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guide



# Arterial Catheterization

