



REF NA-04220-X1A

20<sup>Ga.</sup> 3.81<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guide

# QuickFlash<sup>®</sup> Radial Artery Catheterization

**Contents:**

- 1: QuickFlash<sup>®</sup> Radial Artery Catheter: 20 Ga. x 1-1/2" (3.81 cm) Radiopaque Polyurethane over 21 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: SafetyGlide<sup>™</sup> 25 Ga. x 1" (2.54 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep<sup>®2</sup> One-Step Solution with Hi-Lite Orange<sup>™</sup> Tint
- 1: SharpsAway<sup>®</sup> Disposal Cup
- 1: SharpsAway<sup>®</sup> II Locking Disposal Cup
- 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Checklist Sheet
- 1: Medication Label: 1% Lidocaine
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm<sup>®3</sup> 10 cm x 12 cm
- 1: Dressing: BIOPATCH<sup>®4</sup>
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Tubing<sup>±</sup>: 6" Extension with One-Way Stopcock

<sup>1</sup>A trademark of Becton, Dickinson and Company.  
<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.  
<sup>3</sup>A registered trademark of 3M Company.  
<sup>4</sup>A registered trademark of Johnson & Johnson Corporation.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**

**California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**±Precaution: Contains Phthalate: DEHP.**

**Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.**

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).


Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



LOT



Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex<sup>®</sup>**  
[www.teleflex.com](http://www.teleflex.com)

LBL059319 R01 (2022-02)



REF NA-04220-X1A

20<sup>Ga.</sup> 3.81<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guide

# QuickFlash<sup>®</sup> Radial Artery Catheterization