



REF NA-04550-K

(IPNXXXXXX)

20 Ga.

12 cm catheter length

.025 inch dia. spring-wire guide

Arterial Catheterization

Contents:

- 1: One-Lumen Catheter: 20 Ga. (1.2 mm OD) x 12 cm Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other) Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: Injection Needle: SafetyGlide™¹ 25 Ga. x 1" (2.54 cm) Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: Injection Needle: SafetyGlide™¹ 23 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 2: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep®² One-Step Solution with Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: Drape: 41" x 55" with 4" fenestration Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: Towel: 13" x 18" (33 cm x 46 cm) Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 1: Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) Origin: XX XXXX XX XXX XXXXXXXXXXXXXXXXXXXX

¹A trademark of Becton, Dickinson and Company.
²A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 15 – 30°C (59 – 86°F). Avoid freezing and excessive heat above 40°C (104°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Set Product of xxx
Packaged in xxx

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA



www.teleflex.com

LBL059320 R03 (2024-11)



REF NA-04550-K

20 Ga.

12 cm catheter length

.025 inch dia. spring-wire guide



Arterial Catheterization

