

Dimensions: 4-5/8" x 13-15/16"

PB-10009-000A



ARROW

REF AK-04020-SPC (IPNXXXXXX)

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

### Radial Artery Catheterization with Integral Needle Protection

Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>1</sup> One-Step Frepp® Solution
- 1: SharpsAway® II Locking Disposal Cup
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Towel
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU**

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

K-04020-100C Dimensions: 4 9/16" x 12 21/32 "



ARROW

REF AK-04020-SPC (IPNXXXXXX)

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

### Radial Artery Catheterization with Integral Needle Protection

Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 1.5 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>1</sup> One-Step Frepp® Solution
- 1: SharpsAway® II Locking Disposal Cup
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Towel
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU**

California Prop. 65



⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Product of xxx  
Packaged in xxx

Teleflex®  
www.teleflex.com

LBL059338 R02 (2026-03)

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Product of xxx  
Packaged in xxx

Teleflex®  
www.teleflex.com

LBL059338 R02 (2026-03)



REF AK-04020-SPC

ARROW

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

Radial Artery Catheterization with Integral Needle Protection



REF AK-04020-SPC

ARROW

20 Ga. 4.45 cm catheter length .018 inch dia. spring-wire guide

Radial Artery Catheterization with Integral Needle Protection