



ARROW®

REF AK-04122

22<sup>Ga.</sup> 3.49<sup>cm</sup> catheter length .015<sup>inch dia.</sup> spring-wire guides

## Radial Artery Catheterization

- Contents:**
- 1: Spring-Wire Guide/Tube Assembly: .015" (0.38 mm) dia. x 4-7/8" (12.4 cm)
  - 1: Catheter: 22 Ga. x 1-3/8" (3.49 cm) Radiopaque over 23 Ga. TW Introducer Needle
  - 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
  - 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
  - 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
  - 1: Syringe: 3 mL Luer-Lock
  - 1: Foil Pack 70% Alcohol Prep Pad
  - 1: Skin Protectant Prep
  - 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
  - 1: 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>1</sup> One-Step Frepp® Solution
  - 1: Towel: 13" x 18" (33 cm x 46 cm)
  - 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**  
**California Prop. 65**  
⚠ **WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)  
**Not made with natural rubber latex.**  
Store between 20 – 25°C (68 – 77°F).  
Sterilized using ethylene oxide unless otherwise indicated in the contents list.  
Fluid path components are non-pyrogenic.



ARROW®

REF AK-04122

22<sup>Ga.</sup> 3.49<sup>cm</sup> catheter length .015<sup>inch dia.</sup> spring-wire guides

## Radial Artery Catheterization

- Contents:**
- 1: Spring-Wire Guide/Tube Assembly: .015" (0.38 mm) dia. x 4-7/8" (12.4 cm)
  - 1: Catheter: 22 Ga. x 1-3/8" (3.49 cm) Radiopaque over 23 Ga. TW Introducer Needle
  - 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
  - 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
  - 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
  - 1: Syringe: 3 mL Luer-Lock
  - 1: Foil Pack 70% Alcohol Prep Pad
  - 1: Skin Protectant Prep
  - 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
  - 1: 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>1</sup> One-Step Frepp® Solution
  - 1: Towel: 13" x 18" (33 cm x 46 cm)
  - 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**  
**California Prop. 65**  
⚠ **WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)  
**Not made with natural rubber latex.**  
Store between 20 – 25°C (68 – 77°F).  
Sterilized using ethylene oxide unless otherwise indicated in the contents list.  
Fluid path components are non-pyrogenic.



LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

LBL059343 R00 (2021-11)



LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

LBL059343 R00 (2021-11)



REF AK-04122

ARROW®

22<sup>Ga.</sup> 3.49<sup>cm</sup> catheter length .015<sup>inch dia.</sup> spring-wire guides

## Radial Artery Catheterization



REF AK-04122

ARROW®

22<sup>Ga.</sup> 3.49<sup>cm</sup> catheter length .015<sup>inch dia.</sup> spring-wire guides

## Radial Artery Catheterization