



REF ASK-04020-MIH1

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guide

# Radial Artery Catheterization

**Contents:**

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
  - 1: Injection Needle: SafetyGlide™<sup>1</sup> 25 Ga. x 1" (2.54 cm)
  - 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
  - 1: Syringe: 3 mL Luer-Lock
  - 1: Foil Pack 70% Alcohol Prep Pad
  - 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
  - 2: 1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Frepp® Solution
  - 1: Towel: 13" x 18" (33 cm x 46 cm)
  - 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
  - 1: Needle Holder
  - 1: Scissors
  - 1: Medication Label: 1% Lidocaine
  - 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
  - 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
  - 1: Dressing: OpSite®<sup>3</sup> 2-3/4" x 2-1/4" (7 cm x 6 cm)
  - 1: Dressing: BIOPATCH®<sup>4</sup>
  - 1: Suture: 3-0 Silk with Curved Needle
- <sup>1</sup>A trademark of Becton, Dickinson and Company.
- <sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.
- <sup>3</sup>A registered trademark of T.J. Smith & Nephew, Limited.
- <sup>4</sup>A registered trademark of Johnson & Johnson Corporation.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm.  
www.P65Warnings.ca.gov

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
www.teleflex.com

LBL059356 R01 (2022-02)



REF ASK-04020-MIH1

20<sup>Ga.</sup> 4.45<sup>cm</sup> catheter length .018<sup>inch dia.</sup> spring-wire guide

# Radial Artery Catheterization