



REF ASK-04510-HFH

20 Ga.

12 cm catheter length

.025 inch dia. spring-wire guide

Arterial Line

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow
- Injection Needle: SafetyGlide™1 18 Ga. x 1-1/2" (3.81 cm) Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm)

- Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Syringe: 3 mL Luer-Lock 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- s mt. 1% Lidocaine HCI Solution and Alcohol Prep 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution with Hi-Lite Orange™ Tint SharpsAway® II Locking Disposal Cup SharpsAway® Disposal Cup Maximal Barrier Drape™ with 4" fenestration Needle Holder

- Needle Holder
- Forceps
- 1: Safety Scalpel #11

- Gauze Pad: 2" x 2" (5 cm x 5 cm)
 Gauze Pad: 4" x 4" (10 cm x 10 cm)
 Dressing: Tegaderm®3 I.V. 11.5 cm x 8.9 cm
 Surgical Apparel: Impervious Gown

- Surgical Apparel: Bouffant Cap Surgical Apparel: Mask with Eye Shield Suture: 3-0 Silk with Curved Needle
- ¹A trademark of Becton, Dickinson and Company. ²A registered trademark of CareFusion or one of its subsidiaries.
- ³A registered trademark of 3M Company.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex. Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in

the contents list.

Fluid path components are non-pyrogenic.









Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL059359 R01 (2023-05)









