



REF HF-04020-1-S

20^{Ga.} 3.81^{cm catheter length} .025^{inch dia. spring-wire guide}

Radial Artery Catheterization

Contents:

- 1: Indwelling Catheter: 20 Ga. x 3.81 cm Radiopaque Polyurethane with Integral Suture Wing
 - 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
 - 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
 - 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
 - 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
 - 1: Syringe: 3 mL Luer-Lock
 - 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
 - 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep^{®1} One-Step Solution
 - 1: SharpsAway[®] Disposal Cup
 - 1: SharpsAway[®] II Locking Disposal Cup
 - 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
 - 1: Safety Scalpel: #11
 - 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
 - 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)
 - 1: Suture: 3-0 Braided Silk with Straight Needle
- ¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex.


Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex[®]
www.teleflex.com

LBL059360 R00 (2022-03)



REF HF-04020-1-S

20^{Ga.} 3.81^{cm catheter length} .025^{inch dia. spring-wire guide}

Radial Artery Catheterization