



REF HF-04510-1

20<sup>Ga.</sup>

12<sup>cm catheter length</sup>

.025<sup>inch dia. spring-wire guide</sup>


Arterial Line

Contents:

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>1</sup> One-Step Solution
- 1: SharpsAway® II Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Safety Scalpel: #11
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only  
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
California Prop. 65  
⚠ **WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)  
Not made with natural rubber latex.  
Store between 20 – 25°C (68 – 77°F).  
Sterilized using ethylene oxide unless otherwise indicated in the contents list.  
Fluid path components are non-pyrogenic.



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA



LBL059361 R01 (2022-06)



REF HF-04510-1

20<sup>Ga.</sup>

12<sup>cm catheter length</sup>

.025<sup>inch dia. spring-wire guide</sup>

Arterial Line