

## **Radial Artery Catheterization with Integral Needle Protection**

## Contents

- Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Needle With Integral .018" (0.46 mm) dia. Spring-Wire Guide Injection Needle: SafetyGilde<sup>TM1</sup> 25 Ga. x 1" (2.54 cm) Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Syringe: 3 mL Luer-Lock 5 mL 1% Lidocaine HCI Solution and Alcohol Prep 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution with Hi-Lite Orange<sup>TM</sup> Tint SharpsAway® II Locking Disposal Cup Drape: 36" x 41" (31 cm x 16 cm) Needle Holder Safety Scalpei: #11 Medication Label: 1% Lidocaine Gauze Pad: 2", x 2" (5 cm x 5 cm) 1:

- 1: 1:
- 1:

- Medication Label: 1% Lidocaine
  Gauze Pad: 2" x2" (5 cm x 5 cm)
  Gauze Pad: 4" x 4" (10 cm x 10 cm)
  Dressing: Tegaderm<sup>®3</sup> 10 cm x 12 cm
  Surgical Apparel: Bouffant Cap
  Surgical Apparel: Mask with Eye Shield
  Suture: 3-0 Silk with Curved Needle

## Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU California Prop. 65 & WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov Not made with natural rubber latex. Store between 20 - 25°C (68 - 77°F) Sterilized using ethylene oxide unless otherwise indicated in the contents list. Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

LOT  $\sum$ 

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

**ARROW** 



REF NA-04020-1ASP

20<sup>Ga.</sup> 4.45<sup>cm</sup><sub>catheter</sub> .018<sup>inch dia.</sup>

ARROW

**Radial Artery Catheterization with Integral Needle** Protection