



5 cm catheter length .021 inch dia. spring-wire guide

Arterial Catheterization

Contents

- Catheter: 20 Ga. x 5 cm Radiopaque PEBA Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)

- spinig-white Guide, Markett. 0.21 (0.55 hinh) dia. x 13-3/4 (35 din) (Straigh Soft Tip on Both Ends)
 Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
 Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm)
 Syringe: 3 mL Luer-Lock
 SharpsAway® Il Locking Disposal Cup
 SharpsAway® Il Locking Disposal Cup
 SharpsAway® Il Locking Disposal Cup
 Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
 Towel: 13" x 18" (33 cm x 46 cm)
 Needle Holder
 Safety Scalpel: #11
 Checklist Sheet
 Gauze Pad: 2" x 2" (5 cm x 5 cm)
 Gauze Pad: 4" x 4" (10 cm x 10 cm)
 Dressing: Tegaderm®2 10 cm x 12 cm
 Surgical Apparel: Bouffant Cap
 Surgical Apparel: Mask with Eye Shield
 Suture: 3-0 Silk with Curved Needle
 trademark of Becton, Dickinson and Company.

- ¹A trademark of Becton, Dickinson and Company. ²A registered trademark of 3M Company.

Catheter Priming Volume: 0.17 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
California Prop. 65

\(\triangle \tr











Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL059377 R00 (2022-02)







