



REF **SAC-01220-X1A**

(IPNXXXXXX)

20 Ga.

12 cm catheter length

.021 inch dia. spring-wire guide



Arterial Catheterization

Contents:

- | | | |
|--|---|--|
| <ul style="list-style-type: none"> 1: Catheter: 20 Ga. x 12 cm Radiopaque PEBA Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Introducer Needle: 20 Ga. x 2-3/4" (7 cm) XTW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Injection Needle: SafetyGlide™¹ 25 Ga. x 1" (2.54 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Injection Needle: SafetyGlide™¹ 23 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Syringe: 5 mL Luer-Slip Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®² One-Step Solution with Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: SharpsAway® Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Maximal Barrier Drape™ with 4" fenestration Origin: XX XXXX XX XXX | <ul style="list-style-type: none"> XXXXXX XXXXXXXXX 1: Towel: 13" x 18" (33 cm x 46 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Needle Holder Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Safety Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Checklist Sheet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Medication Label: 1% Lidocaine Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Dressing: Tegaderm®³ 10 cm x 12 cm Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Surgical Apparel: Impervious Gown Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Dressing: BIOPATCH®⁴ Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Surgical Apparel: Bouffant Cap Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Surgical Apparel: Mask with Eye Shield Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX 1: Suture: 3-0 Silk with Curved Needle Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX | <ul style="list-style-type: none"> 1: HemoHopper® Fluid Receptacle Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXX <p>¹A trademark of Becton, Dickinson and Company.
 ²A registered trademark of CareFusion or one of its subsidiaries.
 ³A registered trademark of 3M Company.
 ⁴A registered trademark of Johnson & Johnson Corporation.</p> |
|--|---|--|

Catheter Priming Volume: 0.20 mL Priming volume is approximate.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).
 Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

