

REF **SAC-01618-X1A**

(IPNXXXXXX)

18 Ga.**16** cm catheter length**.025** inch dia. spring-wire guide

ARROW
ErgoPack



Arterial Catheterization

Contents:

- | | |
|--|---|
| 1: Catheter: 18 Ga. x 16 cm Radiopaque PEBA | 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | 1: Dressing: Tegaderm ^{®3} 10 cm x 12 cm |
| 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW | 1: Surgical Apparel: Impervious Gown |
| 1: Injection Needle: SafetyGlide ^{™1} 25 Ga. x 1" (2.54 cm) | 1: Dressing: BIOPATCH ^{®4} |
| 1: Injection Needle: SafetyGlide ^{™1} 22 Ga. x 1-1/2" (3.81 cm) RW | 1: Surgical Apparel: Bouffant Cap |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: Surgical Apparel: Mask with Eye Shield |
| 1: Syringe: 3 mL Luer-Lock | 1: Suture: 3-0 Silk with Curved Needle |
| 1: Syringe: 5 mL Luer-Slip | 1: HemoHopper [®] Fluid Receptacle |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | ¹ A trademark of Becton, Dickinson and Company. |
| 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®2} One-Step Solution with Hi-Lite Orange [™] Tint | ² A registered trademark of CareFusion or one of its subsidiaries. |
| 1: SharpsAway [®] Disposal Cup | ³ A registered trademark of 3M Company. |
| 1: SharpsAway [®] II Locking Disposal Cup | ⁴ A registered trademark of Johnson & Johnson Corporation. |
| 1: Maximal Barrier Drape [™] with 4" fenestration | |
| 1: Towel: 13" x 18" (33 cm x 46 cm) | |
| 1: Needle Holder | |
| 1: Safety Scalpel: #11 | |
| 1: Checklist Sheet | |
| 1: Medication Label: 1% Lidocaine | |
| 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) | |

Catheter Priming Volume: 0.24 mL Priming volume is approximate.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

