

REF **SAC-01220-1A**

(IPNXXXXXX)

20 Ga.**12** cm
catheter
length**.021** inch dia.
spring-wire
guide

Arterial Catheterization



Contents:

- | | |
|--|--|
| <ul style="list-style-type: none"> 1: Catheter: 20 Ga. x 12 cm Radiopaque PEBA 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends) 1: Introducer Needle: 20 Ga. x 2-3/4" (7 cm) XTW 1: Injection Needle: SafetyGlide™¹ 25 Ga. x 1" (2.54 cm) 1: Injection Needle: SafetyGlide™¹ 22 Ga. x 1-1/2" (3.81 cm) RW 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) 1: Syringe: 3 mL Luer-Lock 1: Syringe: 5 mL Luer-Slip 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep 2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®² One-Step Solution with Hi-Lite Orange™ Tint 1: SharpsAway® Disposal Cup 1: SharpsAway® II Locking Disposal Cup 1: Drape: 41" x 55" with 4" fenestration 1: Towel: 13" x 18" (33 cm x 46 cm) 1: Needle Holder 1: Safety Scalpel: #11 1: Medication Label: 1% Lidocaine 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) | <ul style="list-style-type: none"> 1: Dressing: Tegaderm®³ 10 cm x 12 cm 1: Surgical Apparel: Impervious Gown 1: Surgical Apparel: Bouffant Cap 1: Surgical Apparel: Mask with Eye Shield 1: Suture: 3-0 Silk with Curved Needle <p>¹A trademark of Becton, Dickinson and Company.
 ²A registered trademark of CareFusion or one of its subsidiaries.
 ³A registered trademark of 3M Company.</p> |
|--|--|

Catheter Priming Volume: 0.20 mL Priming volume is approximate.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

