



REF SAC-01618-1A

18 Ga. 16 cm catheter length .025 inch dia. spring-wire guide

# Arterial Catheterization

Contents:

- 1: Catheter: 18 Ga. x 16 cm Radiopaque PEBA

1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW

1: Injection Needle: SafetyGlide™<sup>1</sup> 25 Ga. x 1" (2.54 cm)

1: Injection Needle: SafetyGlide™<sup>1</sup> 22 Ga. x 1-1/2" (3.81 cm) RW

1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

1: Syringe: 3 mL Luer-Lock

1: Syringe: 5 mL Luer-Slip

1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep

2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint

1: SharpsAway® Disposal Cup

1: SharpsAway® II Locking Disposal Cup

1: Drape: 41" x 55" with 4" fenestration

1: Towel: 13" x 18" (33 cm x 46 cm)

1: Needle Holder

1: Safety Scalpel: #11

1: Medication Label: 1% Lidocaine

4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®<sup>3</sup> 10 cm x 12 cm

1: Surgical Apparel: Impervious Gown

1: Surgical Apparel: Bouffant Cap

1: Surgical Apparel: Mask with Eye Shield

1: Suture: 3-0 Silk with Curved Needle

<sup>1</sup>A trademark of Becton, Dickinson and Company.

<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

<sup>3</sup>A registered trademark of 3M Company.

Catheter Priming Volume: 0.24 mL Priming volume is approximate.

Not made with natural rubber latex.  
Store between 20 – 25°C (68 – 77°F).  
Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

