





## REF EU-25052-HP

50 cm catheter length .018 inch dia. spring-wire guide

## Pressure Injectable Two-Lumen PICC

## Contents:

- Two Lumen Catheter: 5 Fr. (1.81 mm OD) x 50 cm, Pressure Injectable, T-
- Port Connector, Blue FlexTip®, Contamination Guard and Placement Wire GlideThru™ Peel-Away Sheath: 5 Fr. x 2-3/4" (7 cm) Radiopaque over 5 Fr. Dilator
- Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End Straight Stiff Tip on Other) with Arrow
- Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW Syringe: 10 mL Luer-Lock
- Dust Cap: Non-Vented
- SecondSite™ Adjustable Hub: Fastener SecondSite™ Adjustable Hub: Catheter Clamp
- SharpsAway® II Locking Disposal Cup
- Catheter Trimmer
- Safety Scalpel: #11
- Patient Information Booklet Chart Sticker
- Paper Tape Measure

contraindications exist.

- Dressing: GRIP-LOK®1 Catheter Fixation Device

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Not made with natural rubber latex. Fluid path components are non-pyrogenic.

Contraindications: The Pressure Injectable PICC is contraindicated wherever there is presence of device related infections or presence of Clinical assessment of the patient must be completed to ensure no

2:	Valve: MicroCLAVE™2 Neutral Displacement Connector
1A	registered trademark of TIDI.
2A	trademark of ICU Medical, Inc.

All components are CE 2797 unless otherwise noted.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate†† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (18 Ga.)	0.4	356	3106	4
Proximal (18 Ga.)	0.43	396	3490	4

<sup>††</sup> Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

\*\* Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.





LOT


















Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL059707 R02 (2023-09)









50 cm catheter length .018 inch dia. spring-wire guide



 $<sup>\</sup>dagger$  Flow rate values are approximate and are determined using deionized water at 100 cm head height.