## ARROW



## **REF CS-16123-F**

 $12^{\text{Fr.}} 3^{\text{Lumen}}$ 

25<sup>cm</sup><sub>catheter</sub> .035<sup>inch dia.</sup> spring-wire guide

## **Three-Lumen CVC**

Contents:

- Three-Lumen Catheter: 12 Fr. (4.1 mm OD) x 25 cm Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on 1: 1:
- One End "J" Tip on Other) with Arrow Advancer Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle Pressure Transduction Probe 1:
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm 1:
- 1:
- Dust Cap: Non-Vented Drape: 24" x 36" (61 cm x 91 cm) Scalpel: #11 3: 1:
- 1:


All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (16 Ga.)	0.41	2200	
Medial (12 Ga.)	1.5	9288	
Proximal (12 Ga.)	1.5	10292	

\* Priming volumes are approximate and are done without accessories. † Flow rate values are approximate and are determined using deionized water at 100 cm head height.



en be bg hr cs da nl et fi fr de el hu is it lv lt



ECREP Teleflex Medical IDA Business and Technology Park Dublin Road, Athlone 2797 Co. Westmeath, Ireland

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL060271 R01 (2023-10)

**ARROW** 



LOT

 $\leq$ 

## **Three-Lumen CVC**

[]





