ARROW

REF PL-01055

Arrow[®] GlideThru[®] Peel-Away Sheath over Dilator

5.5^{Fr.} 10^{cm}

Dilator accepts .018" dia. or smaller Spring-Wire Guide. Sheath accepts 5.5 Fr. or smaller catheter.

Contents:

1: GlideThru™ Peel-Away Sheath: 5.5 Fr. x 4" (10 cm) Radiopaque over 5.5 Fr. Dilator

ARROW

REF PL-01055

Arrow[®] GlideThru[®] Peel-Away Sheath over Dilator



Dilator accepts .018" dia. or smaller Spring-Wire Guide. Sheath accepts 5.5 Fr. or smaller catheter.

Contents:

1: GlideThru™ Peel-Away Sheath: 5.5 Fr. x 4" (10 cm) Radiopaque over 5.5 Fr. Dilator

ARROW

REF PL-01055

Arrow[®] GlideThru[®] Peel-Away Sheath over Dilator



Dilator accepts .018" dia. or smaller Spring-Wire Guide. Sheath accepts 5.5 Fr. or smaller catheter.

Contents:

1: GlideThru™ Peel-Away Sheath: 5.5 Fr. x 4" (10 cm) Radiopaque over 5.5 Fr. Dilator

Fluid path components are non-pyrogenic.







Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL061003 R00 (2022-07)