

DO NOT BEND

ARROW®



REF SAC-01620 (IPNXXXXXX)

20^{Ga.}16^{cm catheter length}

Arterial Catheterization Set

Contents:

- 1: Indwelling Catheter: 20 Ga. x 16 cm Radiopaque Polyurethane
 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 19-3/4" (50 cm) (Straight Soft Tip on Both Ends) with Arrow Advancer
 1: Introducer Needle: 20 Ga. x 2-3/4" (7 cm) XTW

Catheter Priming Volume: 0.21 mL Priming Volume is approximate.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

Contraindications: Inadequate collateral arterial circulation.

- cs** Arteriální katetrizační set
da Arterielt kateteriseringssæt
nl Arteriële katheterisatieset
et Arterikanüüli komplekt
fi Valtimokatetrisarja
fr Jeu de cathétérisme artériel
de Besteck zur Katheterisierung von Arterien
el Σετ αρτηριακού καθετηριασμού
hu Artériás katéter készlet
it Set per cateterismo arterioso
lt Arterijų kateterizavimo rinkinys
no Arteriell kateteriseringssett
pl Zestaw do cewnikowania tętnic
pt Conjunto de cateterização arterial
ro Set de cateterizare arterială
ru Набор для артериальной катетеризации
sk Arteriálna katetrizačná sada
sl Pribor za arterijsko kateterizacijo
es Conjunto de cateterización arterial
sv Artärkateteriseringsset
tr Arteriyel Kateterizasyon Seti



(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (10)LotNumber



LotNumber



YYYY-MM-DD

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL062141 (2021-06)



Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland



REF SAC-01620

Arterial Catheterization Set

20^{Ga.}16^{cm catheter length}

DO NOT BEND

ARROW®