ARRC









(IPNXXXXXX)

6 cm catheter length

Pressure Injectable Arrowg+ard Blue Plus® Four-Lumen CVC

Contents:

- Arrowg+ard Blue Plus®1 Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 16 cm, Pressure Injectable
- Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End "J" Tip on Other) with Arrow Advancer with ECG Mark Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- - Pressure Transduction Probe
- Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL Tissue Dilator: 9 Fr. x 10.2 cm
- Dust Cap: Non-Vented
 SecondSite™ Adjustable Hub: Fastener
 SecondSite™ Adjustable Hub: Catheter Clamp
- SharpsAway® II Locking Disposal Cup
- Safety Scalpel: #11

¹Licensed under US Patent No. 7,329,412

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.







Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (16 Ga.)	0.38	2781	5
Medial 1 (14 Ga.)	0.45	4980	10
Medial 2 (18 Ga.)	0.35	1453	not rated
Proximal (18 Ga.)	0.38	1619	not rated

* Priming volumes are approximate and are done without accessories



All components are CE 2797 unless otherwise noted.

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA


Product of xxx Packaged in xxx


(01)XXXXXXXXXXXXXXX (11)YYMMDD


(10)Lot Number

)feleflex® www.teleflex.com

LBL062490 R02 (2025-05)
















REF EU-42854-HPS

Lumen

6 cm catheter length

.032 inch dia. spring-wire guide





[†] Flow rate values are approximate and are determined using deionized water at 100 cm head

Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 152 cm administration set tubing.