



REF ASK-09803-SRH

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

# PSI Kit for use with 7 Fr. Catheters

**Contents:**

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide™<sup>1</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- 1: Injection Needle: SafetyGlide™<sup>1</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 1: Dust Cap: Non-Vented
- 1: SharpsAway® II Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Maximal Barrier Drape™ with 4" fenestration
- 1: Drape: 29" x 18" (74 cm x 46 cm) with adhesive
- 1: Safety Scalpel: #11
- 1: Stopcock: 4-way High-Flow
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>A trademark of Becton, Dickinson and Company.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**Caution: This Product Contains Dry Natural Rubber.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
 www.teleflex.com

LBL063872 R00 (2022-02)



REF ASK-09803-SRH

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

# PSI Kit for use with 7 Fr. Catheters