

Dimensions: 6-3/8" (16 cm) x 8-1/4" (21 cm)



REF **PI-01351-LS** (IPNXXXXXX)

for use with Arrow

3 Fr. PICC

7 cm

.016 inch dia. spring-wire guide

Insertion Set for use with Arrow® 3 Fr. PICC

Contents:

- 1: GlideThru™ Peel-Away Sheath: 3 Fr. x 2-3/4" (7 cm) Radiopaque over 3 Fr. Dilator
origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Spring Wire Guide Assembly: .016" (0.40 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip On Both Ends) origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW origin:xxxxxxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW origin:xxxxxxxxxxxxxxxxxxxxxxxx

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.



CAS # 7440-48-4

EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Set Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL064451 R03 (2024-05)