



REF AK-04510-S

20<sup>Ga.</sup>

12<sup>cm</sup>  
catheter  
length

.025<sup>inch dia.</sup>  
spring-wire  
guide

# Arterial Line

**Contents:**

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 2: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep®<sup>1</sup> One-Step Solution
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Scalpel: #11
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**

**California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.**

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

**Due To A Temporary Supply Shortage  
A Clear ChloroPrep Has Been Added To This Kit To Replace  
The Orange Tinted ChloroPrep.**



LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL067091 R00 (2022-06)



REF AK-04510-S

20<sup>Ga.</sup>

12<sup>cm</sup>  
catheter  
length

.025<sup>inch dia.</sup>  
spring-wire  
guide

# Arterial Line