



Arterial Catheterization

Contents:

- Contents:

 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide

 1: Injection Needle: SafetyGlide **1 25 Ga. x 1" (2.54 cm)

 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

 1: Syringe: 3 mL Luer-Lock

 1: 5 mL 1% Lidocaine HCI Solution and Alcohol Prep

 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution

 1: Wing Clip: Suture

 1: SharpsAway® Disposal Cup

 1: SharpsAway® IL Locking Disposal Cup

 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration

 1: Towel: 13" x 18" (33 cm x 46 cm)

 1: Needle Holder

 1: Safety Scalpel: #11

 1: Medication Label: 1% Lidocaine

 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)

 1: Dressing: Tegaderm®3 10 cm x 12 cm

 1: Surgical Apparel: Bouffant Cap

 1: Surgical Apparel: Mask with Eye Shield

- 1: Suture: 3-0 Silk with Curved Needle
- 1. Souther 3.0 oils with Couver Needle 2.1. A trademark of Becton, Dickinson and Company.

 2A registered trademark of CareFusion or one of its subsidiaries.

 3A registered trademark of 3M Company.

Store between 20 – 25°C (68 – 77°F).
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
Fluid path components are non-pyrogenic.
Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Due To A Temporary Supply Shortage A Clear ChloraPrep Has Been Added To This Kit To Replace The Orange Tinted ChloraPrep.

























Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA



LBL067092 R00 (2022-06)







