16 cm catheter length .025 inch dia. spring-wire quide

















Scalpel: #11

REF SA-04018

Femoral Artery Catheterization

Femoral Artery Catheterization

injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.

Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)

1: Indwelling Catheter: 18 Ga. x 16 cm with Integral Suture Wing

Indwelling Catheter: 18 Ga. x 16 cm with Integral Suture Wing

Introducer Needle: 18 Ga. x 3" (7.62 cm) RW

1: Suture: 2-0 Silk with Curved Needle

- Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- Introducer Needle: 18 Ga. x 3" (7.62 cm) RW
- 1: Scalpel: #11
- 1: Suture: 2-0 Silk with Curved Needle

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient

16 cm catheter length .025 inch dia. spring-wire quide





Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA











Femoral Artery Catheterization

Contents:

- 1: Indwelling Catheter: 18 Ga. x 16 cm with Integral Suture Wing
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 18 Ga. x 3" (7.62 cm) RW

16 cm catheter length .025 spring-wire quide



Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA











1: Suture: 2-0 Silk with Curved Needle

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.