

ARROW®



REF SS-14702

7 Fr. 16 cm

**Two-Lumen CVC for use only with Arrow Percutaneous Sheath Introducer (8.5 Fr. and 9 Fr.)**

Contents:

- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) Catheter Contamination Guard
- 2: Dust Cap: Non-Vented

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**  
[www.teleflex.com/IFU](http://www.teleflex.com/IFU)

**Not made with natural rubber latex.**  
Fluid path components are non-pyrogenic.



Lumen

Priming Volume\* (mL)

Gravity Flow Rate† (mL/hr)

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.35	3287
Proximal (16 Ga.)	0.39	3453

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL067454 R01 (2022-09)

**Two-Lumen CVC for use only with Arrow Percutaneous Sheath Introducer (8.5 Fr. and 9 Fr.)**



REF SS-14702

7 Fr. 16 cm

ARROW®