

ARROW®



REF SS-14703

7 Fr. 16 cm

Three-Lumen CVC for use only with Arrow Percutaneous Sheath Introducer (8.5 Fr. and 9 Fr.)

Contents:

- 1: Three-Lumen Catheter: 7 Fr. (2.5 mm OD) Contamination Guard (14 cm Catheter Length when used with Percutaneous Sheath Introducer)
- 3: Dust Cap: Non-Vented

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen

Priming Volume* (mL)

Gravity Flow Rate† (mL/hr)

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.34	2913
Medial (18 Ga.)	0.27	1403
Proximal (18 Ga.)	0.27	1502

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



LOT



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL067468 R01 (2022-09)

Three-Lumen CVC for use only with Arrow Percutaneous Sheath Introducer (8.5 Fr. and 9 Fr.)



REF SS-14703

7 Fr. 16 cm

ARROW®