
PB-10000-850 8.5 X 10.5



- SecondSite™ Adjustable Hub: Catheter Clamp 1.
- Catheter Trimmer
- Safety Scalpel: #11
- Patient ID Card 1.
- Chart Sticker
- Patient Information Booklet
- Checklist/CLIP Sheet
- Paper Tape Measure 1:
- Dressing: STATLOCK®1 PICC Plus Catheter Stabilization Device 1:
- ¹A registered trademark of C. R. Bard, Inc.
- All components are CE 2797 unless otherwise noted.

O Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate†† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (18 Ga.)	0.52	416	3660	5
Proximal (18 Ga.)	0.53	409	3560	5

† Flow rate values are approximate and are determined using deionized water at 100 cm head

** Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum

++ Pump flow rates are determined at maximum pump pressure of 10 psig and represent

* Priming volumes are approximate and are done without accessories.

using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.

height.

approximate flow capabilities.

Not made with natural rubber latex

Contraindications: The Pressure Injectable PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.

//CAS # 7440-48-4

LBL068049 R01 (2025-01)

