





(IPNXXXXXX)

2 Lumen 11.5 cm device length

.035 inch dia. spring-wire quide

## MAC™ Two-Lumen Central Venous Access for use with 7.5 - 8 Fr. Catheters

## Contents:

- Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter Origin: XX XXXX XX XXX XXXXXXXX

- <sup>1</sup>A trademark of Becton, Dickinson and Company. <sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience.

Please ensure that these labels are applied to the correct syringe and corresponding medication.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (9 Fr.) no catheter	1.8	31440
Distal (9 Fr.) with 7 Fr. catheter	1.45	16280
Distal (9 Fr.) with 8 Fr. catheter	1.15	9040
Proximal (12 Ga.) no catheter	0.6	8970

<sup>\*</sup> Priming volumes are approximate and are done without accessories.



Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

**Teleflex**®

LBL069457 R02 (2024-11)



**LOT** Lot Number

YYYY-MM-DD



YYYY-MM-DD

2 Lumen 11.5 cm device length

1035 inch dia. spring-wire





(01)XXXXXXXXXXXXXXX (17)YYMMDD (11)YYMMDD

(10)Lot Number





<sup>†</sup> Flow rate values are approximate and are determined using deionized water at 100 cm head height.