



REF ASK-04510-UPM

20 Ga. 12 cm catheter length .025 inch dia. spring-wire guide

Arterial Catheterization

Contents:

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP

1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other)

1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW

1: Injection Needle: SafetyGlide™<sup>1</sup> 25 Ga. x 1" (2.54 cm)

1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

2: Syringe: 3 mL Luer-Lock

1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep

1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint

1: SharpsAway® II Locking Disposal Cup

1: SharpsAway® Disposal Cup

1: Drape: 41" x 55" with 4" fenestration

1: Needle Holder

1: Scissors

1: Forceps

1: Safety Scalpel: #11

3: Gauze Pad: 2" x 2" (5 cm x 5 cm)

2: Suture: 2-0 Silk with Curved Needle

1: Tubing±: 20" (50 cm) Extension
- <sup>1</sup>A trademark of Becton, Dickinson and Company.

<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

**Rx only**  
**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
**California Prop. 65**  
**⚠ WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)  
**±Precaution:** Contains Phthalate: DEHP.  
**Warning:** The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.  
**Not made with natural rubber latex.**  
Store between 20 – 25°C (68 – 77°F).  
Sterilized using ethylene oxide unless otherwise indicated in the contents list.  
Fluid path components are non-pyrogenic.



LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

LBL069460 R00 (2022-12)



REF ASK-04510-UPM

20 Ga. 12 cm catheter length .025 inch dia. spring-wire guide

Arterial Catheterization