



REF ASK-04510-UPM (IPNXXXXXX)

20 Ga. 12 cm catheter length .025 inch dia. spring-wire guide

Arterial Catheterization

Contents:

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP origin:xxxxxxxxxxxxxxxxxxxx
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW origin:xxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: SafetyGlide™¹ 25 Ga. x 1" (2.54 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 2: Syringe: 3 mL Luer-Lock origin:xxxxxxxxxxxxxxxxxxxx
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep origin:xxxxxxxxxxxxxxxxxxxx
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®² One-Step Solution with Hi-Lite Orange™ Tint origin:xxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® II Locking Disposal Cup origin:xxxxxxxxxxxxxxxxxxxx
- 1: SharpsAway® Disposal Cup origin:xxxxxxxxxxxxxxxxxxxx
- 1: Drape: 41" x 55" with 4" fenestration origin:xxxxxxxxxxxxxxxxxxxx
- 1: Needle Holder origin:xxxxxxxxxxxxxxxxxxxx
- 1: Scissors origin:xxxxxxxxxxxxxxxxxxxx
- 1: Forceps origin:xxxxxxxxxxxxxxxxxxxx
- 1: Safety Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxx
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 2: Suture: 2-0 Silk with Curved Needle origin:xxxxxxxxxxxxxxxxxxxx
- 1: Tubing±: 20" (50 cm) Extension origin:xxxxxxxxxxxxxxxxxxxx

¹A trademark of Becton, Dickinson and Company.
²A registered trademark of CareFusion or one of its subsidiaries.

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
California Prop. 65
⚠ WARNING: Cancer and Reproductive Harm. www.P65Warnings.ca.gov
±Precaution: Contains Phthalate: DEHP.
Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.
Not made with natural rubber latex.
Store between 20 – 25°C (68 – 77°F).
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
Fluid path components are non-pyrogenic.

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

(01)XXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Set Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com
LBL069460 R01 (2024-06)



REF ASK-04510-UPM

20 Ga. 12 cm catheter length .025 inch dia. spring-wire guide



Arterial Catheterization