



CDC-45703-XPN1A

(IPNXXXXXX)

7 Fr.

3 Lumen

20 cm catheter length

.032 inch dia. spring-wire guide

This Product Does Not Contain a 5 mL Arrow® Raulerson Syringe or Pressure Transduction Probe Due to a Supply Shortage, and Contains An Additional 5 mL Syringe.

ARROW®

ErgoPack Complete

Pressure Injectable Arrowg+ard Blue Plus® Three-Lumen CVC



Contents:

- 1: Arrowg+ard Blue Plus®¹ Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm, Pressure Injectable
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: SafetyGlide™² 25 Ga. x 1" (2.54 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Luer-Slip Syringe
- 1: Injection Needle: SafetyGlide™² 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: 10 mL Luer-Lock Pre-Filled Saline Syringe
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®³ One-Step Solution with HI-Lite Orange™ Tint
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® II Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Maximal Barrier Drape™ with 4" fenestration
- 1: Towel
- 1: Needle Holder
- 1: Safety Scalpel: #11

- 1: Checklist/CLIP Sheet
 - 1: Sterile Procedure Sign
 - 1: Medication Label: 1% Lidocaine
 - 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
 - 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
 - 1: Surgical Apparel: Impervious Gown
 - 1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm
 - 1: Dressing: BIOPATCH®⁴
 - 1: Transducer Cover: 14 cm x 147 cm
 - 1: Surgical Apparel: Mask with Eye Shield
 - 1: Surgical Apparel: Bouffant Cap
 - 1: Suture: 3-0 Silk with Curved Needle
 - 1: Tubing±: 20" (50 cm) Extension
 - 1: HemoHopper® Fluid Receptacle
- ¹Licensed under US Patent No. 7,329,412.
²A trademark of Becton, Dickinson and Company.
³A registered trademark of CareFusion or one of its subsidiaries.
⁴A registered trademark of Johnson & Johnson Corporation.

±Precaution: Contains Phthalate: DEHP.

Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



Lumen

Priming
Volume*
(mL)

Gravity
Flow
Rate†
(mL/hr)

MAX Pressure
Injection
Flow Rate**
(mL/sec)

Distal (16 Ga.)	0.39	2283	10
Medial (18 Ga.)	0.39	913	5
Proximal (18 Ga.)	0.40	996	5

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

** Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 152 cm administration set tubing.

