





Arrowg+ard Blue Plus® Two-Lumen CVC

Contents:

Arrowg+ard Blue Plus®¹ Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End-"J" Tip on Other) with Arrow Advancer with ECG Mark Introducer_Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW

Pressure Transduction Probe

Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL

Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm

2.

Dust Cap: Non-Vented SecondSite™ Adjustable Hub: Fastener 1.

SecondSite™ Adjustable Hub: Catheter Clamp 1Licensed under US Patent No. 7,329,412.

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

de el hu is
it Iv It no
pl pt ro ru
sr sk sl es
sv tr uk









6 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	
Distal (14 Ga.)	0.57	4731	_
Proximal (18 Ga.)	0.40	913	

riming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.





EU Authorized Representative and Importer:



Teleflex Medical IDA Business and Technology Park Dublin Road, Athlone 2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

> Product of xxx Packaged in xxx




www.teleflex.com



Lot Number



Y-MM-DD





(01)XXXXXXXXXXXXXXX (17)YYMMDD (11)YYMMDD (10)Lot Number

LBL074106 R01 (2024-09)





Lumen

catheter lenath

32 spring-wire guide





