

DO NOT BEND

ARROW®



REF SAC-01220-PBX
(IPNXXXXXX)

20^{Ga.} 12^{cm catheter length}

Arterial Catheterization Set

Contents:

- 1: Catheter: 20 Ga. x 12 cm Radiopaque PEBA Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Introducer Needle: 20 Ga. x 2-3/4" (7 cm) XTW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.20 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU
California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

en bg cs da nl et fi

fr de el hu it lv lt no

pl pt ro ru sk sl es sv

tr

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
2797

All components are CE 2797 unless otherwise noted.

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Set Product of xxx
Packaged in xxx



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Teleflex®
www.teleflex.com

LBL075266 R00 (2025-05)



REF SAC-01220-PBX

20^{Ga.} 12^{cm catheter length}

Arterial Catheterization Set

ARROW®

DO NOT BEND