



REF CDC-22122-CN1AU (IPNXXXXXX)

12^{Fr.} 2^{Lumen} 16^{cm catheter length} .035^{inch dia. spring-wire guide}

ARROW®

ErgoPack[®]
Complete

Arrowg+ard Blue® You-Bend™ Two-Lumen Hemodialysis Catheterization for High Volume Infusions

UDI



Contents:


- | | |
|---|--|
| 1: Two-Lumen Catheter†: 12 Fr. (4.06 mm OD) x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment, Extension Line Clamps | 1: Safety Scalpel: #11 |
| 1: Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer | 1: Checklist/CLIP Sheet |
| 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle | 1: Sterile Procedure Sign |
| 1: Injection Needle: SafetyGlide™ ¹ 25 Ga. x 1" (2.54 cm) | 1: Medication Label: 1% Lidocaine |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Pressure Transduction Probe | 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe | 1: Surgical Apparel: Impervious Gown |
| 1: Injection Needle: SafetyGlide™ ¹ 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe | 1: Dressing: Tegaderm® ³ CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm) |
| 1: Syringe: 3 mL Luer-Lock | 1: Transducer Cover: 14 cm x 147 cm |
| 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm | 1: Surgical Apparel: Mask with Eye Shield |
| 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm | 1: Surgical Apparel: Bouffant Cap |
| 2: 10 mL Luer-Lock Pre-Filled Saline Syringe | 1: Suture: 3-0 Silk with Curved Needle |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | 1: Tubing‡: 20" (50 cm) Extension |
| 2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep® ² One-Step Solution with Hi-Lite Orange™ Tint | 1: HemoHopper® Fluid Receptacle |
| 1: SharpsAway® II Locking Disposal Cup | |
| 1: SharpsAway® Disposal Cup | |
| 1: Maximal Barrier Drape™ with 4" fenestration | |
| 1: Towel | |
| 1: Needle Holder | |

‡Precaution: Contains Phthalate: DEHP.
Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

‡MR Conditional
Not made with natural rubber latex.
Store between 20 – 25°C (68 – 77°F).


Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels.
The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs.



Priming Volume*	
Distal	1.70 mL
Proximal	1.60 mL

*Priming volumes are approximate and are done without accessories.



† In vitro tests performed using 3:1 glycerine/saline at 37°C.