



REF **DLX-35541-CURVC** (IPNXXXXXX)

4 Fr. **1** Lumen **55** cm catheter length **.018** inch dia. spring-wire guide

ARROW®

Pressure Injectable One-Lumen PICC pre-loaded with Arrow® NaviCurve™ Stylet

UDI




Contents:

- | | |
|--|---|
| 1: One Lumen TaperFree® Catheter,†††: 4 Fr. (1.45 mm OD) x 55 cm, Pressure Injectable, Arrow® NaviCurve™ Stylet and T-Port Connector | 1: Patient Information Booklet |
| 1: GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator | 1: Medication Label: 1% Lidocaine |
| 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer | 2: Paper Tape Measure |
| 1: Injection Needle: Eclipse™1 25 Ga. x 1" (2.54 cm) | 1: Tourniquet |
| 1: Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW | 1: Foam Electrodes (3 per pack) |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 10: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW | 10: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: Syringe: 5 mL Luer-Lock | 1: Surgical Apparel: Impervious Gown |
| 4: 10 mL Luer-Lock Pre-Filled Saline Syringe | 1: Dressing: Tegaderm®3 CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm) |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | 1: Transducer Cover with Gel and Rubber Bands |
| 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution | 1: Remote Cover |
| 2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution with Hi-Lite Orange™ Tint | 1: Dressing: STATLOCK®4 PICC Plus Catheter Stabilization Device |
| 1: Dust Cap: Non-Vented | 1: Surgical Apparel: Bouffant Cap |
| 1: SharpsAway® II Locking Disposal Cup | 1: Surgical Apparel: Mask |
| 1: SharpsAway® Disposal Cup | 1: Surgical Apparel: Mask with Eye Shield |
| 1: Catheter Trimmer | 1: Tape: Steri-Strip®3 |
| 1: Towel: 24" x 36" (61 cm x 91 cm) | 1: T-Piece Cover |
| 1: Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations | |
| 1: Safety Scalpel: #11 | |
| 1: Sterile Procedure Sign | |
| 1: Patient ID Card | |

***MR Safe
Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).
Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The Pressure Injectable PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.



 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate†† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Pink (17 Ga.)	0.63	970	5300	5

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

†† Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

** Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.