



REF **DLX-35541-VPSC** (IPNXXXXXX)

4 Fr. **1** Lumen **55** cm catheter length **.018** inch dia. spring-wire guide

Pressure Injectable One-Lumen PICC pre-loaded with Arrow® VPS Precision™ Stylet

UDI



Contents:

- | | |
|--|--|
| 1: One Lumen TaperFree® Catheter, 4 Fr. (1.45 mm OD) x 55 cm, Pressure Injectable, Contamination Guard, Arrow® VPS Precision™ Stylet and Tuohy-Borst Adapter | 1: Sterile Procedure Sign |
| 1: GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator | 1: Patient ID Card |
| 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer | 1: Patient Information Booklet |
| 1: Stiffening Stylet (0.36mm dia. x 75cm with straight distal tip) | 1: Medication Label: 1% Lidocaine |
| 1: Injection Needle: Eclipse™ 1 25 Ga. x 1" (2.54 cm) | 2: Paper Tape Measure |
| 1: Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW | 1: Tourniquet |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) | 1: Foam Electrodes (3 per pack) |
| 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW | 10: Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: Syringe: 5 mL Luer-Lock | 10: Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 4: 10 mL Luer-Lock Pre-Filled Saline Syringe | 1: Surgical Apparel: Impervious Gown |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep | 1: Dressing: Tegaderm® CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm) |
| 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution | 1: Transducer Cover with Gel and Rubber Bands |
| 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution with Hi-Lite Orange™ Tint | 1: Dressing: STATLOCK® PICC Plus Catheter Stabilization Device |
| 1: Dust Cap: Non-Vented | 1: Surgical Apparel: Bouffant Cap |
| 1: SharpsAway® II Locking Disposal Cup | 1: Surgical Apparel: Mask |
| 1: SharpsAway® Disposal Cup | 1: Surgical Apparel: Mask with Eye Shield |
| 1: Catheter Trimmer | 1: Tape: Steri-Strip® |
| 1: Towel: 24" x 36" (61 cm x 91 cm) | |
| 1: Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations | |
| 1: Safety Scalpel: #11 | |


***MR Safe

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The Pressure Injectable PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.



Lumen

Priming
Volume*
(mL)

Gravity
Flow
Rate†
(mL/hr)

Pump
Flow Rate††
(mL/hr)

MAX Pressure
Injection
Flow Rate**
(mL/sec)

Pink (17 Ga.)	0.63	970	5300	5
---------------	------	-----	------	---

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

†† Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

** Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.