



REF PR-05541-HPX (IPNXXXXXX)

4 Fr. **1** Lumen **55** cm catheter length **.018** inch dia. spring-wire guide

Pressure Injectable One-Lumen PICC with 130 cm Wire

Contents:

- | | |
|--|--|
| 1: One Lumen Catheter: 4 Fr. (1.45 mm OD) x 55 cm, Pressure Injectable, Blue FlexTip® Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Towel: 24" x 36" (61 cm x 91 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: GlideThru™ Peel-Away Sheath: 4 Fr. x 4" (10 cm) Radiopaque over 4 Fr. Dilator Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Safety Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Patient ID Card Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 51-3/16" (130 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Chart Sticker Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Patient Information Booklet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: Injection Needle: Eclipse™ 21 Ga. x 1" (2.54 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Medication Label: 1% Lidocaine Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: Injection Needle: SafetyGlide™ 21 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Medication Label: 0.9% Saline Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Checklist/CLIP Sheet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 2: Paper Tape Measure Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 2: Syringe: 10 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Tourniquet Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 5: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: 10 mL Luer-Lock Pre-Filled Saline Syringe Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 10: Gauze Pad: 4" x 4" (10 cm x 10 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Dressing: I.V. Advanced 8.5 cm x 11.5 cm Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution with Hi-Lite Orange™ Tint Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Dressing: BIOPATCH®³ Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: Dust Cap: Non-Vented Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Dressing: STATLOCK®⁴ PICC Plus Catheter Stabilization Device Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: SecondSite™ Adjustable Hub: Catheter Clamp Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | 1: Tape: Steri-Strip®⁵ Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX |
| 1: SecondSite™ Adjustable Hub: Fastener Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | |
| 1: SharpsAway® II Locking Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | |
| 1: SharpsAway® Disposal Cup Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | |
| 1: Catheter Trimmer Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX | |
- ¹A trademark of Becton, Dickinson and Company.

²A registered trademark of CareFusion or one of its subsidiaries.

³A registered trademark of Johnson & Johnson Corporation.

⁴A registered trademark of C. R. Bard, Inc.

⁵A registered trademark of 3M Company.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

Contraindications: The Pressure Injectable PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.

 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate†† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (18 Ga.)	0.48	436	4075	4

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

†† Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

** Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.



LOT Lot Number

YYYY-MM-DD

 YYYY-MM-DD



(01)XXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Set Product of xxx
Packaged in xxx

 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL080207 R00 (2025-08)



REF PR-05541-HPX

4 Fr. **1** Lumen **55** cm catheter length **.018** inch dia. spring-wire guide



Pressure Injectable One-Lumen PICC with 130 cm Wire